




Publication No. IP 3439
04/96

ILLUSTRATED PARTS LIST

for the

KODAK AUTOMIXER II Plus



PLEASE NOTE

The information contained herein is based on the experience and knowledge relating to the subject matter gained by Kodak prior to publication.

No patent licence is granted by this information

Kodak reserves the right to change this information without notice, and makes no warranty express or implied, with respect to this information. Kodak shall not be liable for any loss or damage, including consequential or special damages, resulting from the use of this information, even if loss or damage is caused by Kodak's negligence or other fault.

TABLE OF CONTENTS

FIGURE:**PAGE**

1. TANKS, LOWER BODIES, COVERS, PANELS, SHROUD, REED CONTACT ARRAYS	
FAN ASSY, CASTERS, BRACKETS	2
2. BOTTOM PLATE (TANKS)	4
3. BASINS	6
4. BASINS (TUBING)	8
5. TUBING	10
6. RINSING	12
7. WATER INLET, DRAIN	14
8. CLAMPS	16
9. ELECTRONIC BOX (TRANSFORMER)	18
10. ELECTRONIC BOX (POWER and REMOTE CONNECTOR)	20
11. ELECTRONIC BOX (SERVICE KEYBOARD, HARNESS, CABLES)	22
12. ELECTRONIC BOX (PRINTED CIRCUIT BOARDS)	24
13. OPERATOR DISPLAY	26
14. BOARDS, CONNECTORS, CABLES	28

APPENDIX:

NUMERICAL INDEX	N-1
ALPHABETICAL INDEX	A-1


**TANKS, LOWER BODIES, COVERS, PANELS, SHROUD, REED CONTACT ARRAYS,
FAN ASSY, CASTERS, BRACKETS**

figure 1

Parts Listing figure 1

Item	Part No.	Description	Quantity
1	927 0106	Tank Developer	1
2	927 0107	Tank Fixer	1
3	927 0136	Lower Body Developer	1
4	927 0126	Lower Body Fixer	1
5	917 1101	Support	2
6	927 0176	Cover	2
7	927 0191	Shaft	4
8	927 0156	Front Panel Assy with Magnet and Screw	1
	917 1451	Magnet	1
	430 2101	Screw	1
9	927 0166	Rear Panel	1
10	927 0116	Shroud	1
11	917 2471	Screw (Phillips)	4
12	448 0271	Washer	4
13	927 1856	Reed Contact Array with Float, Hexagon Nuts, Cable and Connector	2
14	927 1136	Fan Assy including Connector	1
15	430 2063	Screw (Phillips Head, self tapping)	4
16	448 0032	Washer	4
17	917 1381	Cover (Transfer Valves)	4
18	917 0421	Locking Caster or	2
	917 0431	Caster	2
19	448 0181	Washer	4
20	485 2771	Hexagon Nut	4
21	917 6141	Cover (Circulation Pumps)	1
22	917 0681	Valve Bracket	1
23	917 0704	Bracket	1
24	428 1461	Screw	1
25	430 2081	Screw	2

BOTTOM PLATE (TANKS)




figure 2

Parts Listing figure 2

Item	Part No.	Description	Quantity
1	927 0052	Exhaust Tube Assy without O-Ring	2
2	927 0082	Nut	4
3	917 0171	O-Ring	6
4	927 0072	Bottom Plate	2
5	927 0096	Gasket	2
6	485 3481	Hexagon Nut M6	4
7	450 0681	Lock Washer	4
8	927 0061	Bracket	4
9	428 3961	Hexagon Screw M6X40	4

BASINS




figure 3

Parts Listing figure 3

Item	Part No.	Description	Quantity
1	927 0406	Mixing Basin Developer	1
	927 0416	Mixing Basin Fixer	1
CAT	715 7068	Developer + Fixer Basin Assy RP Chemistry (litres)	1
CAT	715 7076	Developer + Fixer Basin Assy RP Chemistry (Gallons)	1
CAT	715 0609	Developer + Fixer Basin Assy RA Chemistry (Gallons)	1
CAT	??? ???? ?	GISD Chemistry (not available yet)	
2	927 0466	Developer Basin without Accessories	1
	927 0456	Fixer Basin without Accessories	1
3	927 0716	Template RP Developer (litres)	1
	927 0706	Template RP Fixer (litres)	1
	927 0736	Template RP Developer (Gallons)	1
	927 0726	Template RP Fixer (Gallons)	1
	927 0756	Template RA Developer (Gallons)	1
	927 0746	Template RA Fixer (Gallons)	1
4	917 0041	Rivet Plastic	8
5	927 0506	Piercer Assy complete with Hexagon Plastic Screw	5
		PN 927 0551 and Tube Nozzle PN 927 0561	
6	927 0282	Coupling Rinse	2
7	927 0281	Screw Cap Rinse	2

BASINS (TUBING)





figure 4

Parts Listing figure 4

Item	Part No.	Description	Quantity
1	927 0551	Hexagon Plastic Screw	5
2	927 0561	Tube Nozzle	5
3	927 0581	T-Fitting	3
4	927 1701	Tube*	4
5	927 1721	Tube*	1
6	927 1711	Tube*	1
7	927 0331	Tube*	2

*** To ensure a constant rinse quality, use original Tubes only !**

TUBING




figure 5

Parts Listing figure 5

Item	Part No.	Description	Quantity
1	927 0041	T - Fitting	2
2	917 0276	Tube Developer (Mixing Tank > Solenoid)	1
3	917 0306	Tube Fixer (Mixing Tank > Solenoid)	1
4	917 6206	Tube Developer/Fixer (Holding Tank > Output)	4
5	927 1736	Tube Fixer (Solenoid > Circulation Pump)	1
6	927 1756	Tube Fixer (Circulation Pump > Mixing Valve)	1
7	917 0266	Tube Developer/Fixer (Mixing Valve > Floating Lid)	2
8	917 1136	Floating Lid	2
9	917 0861	Elbow	2
10	917 0171	O - Ring	4
11	917 0824	Elbow Fitting	2
12	917 0624	Transfer Valve Developer	1
13	917 0634	Transfer Valve Fixer	1
14	927 0496	Circulation Pump	2
15	917 0654	Developer Mixing Valve Assy	1
16	917 0664	Fixer Mixing Valve Assy	1
17	927 1746	Tube Developer (Solenoid > Circulation Pump)	1
18	927 1766	Tube Developer (Circulation Pump > Mixing Valve)	1
19	927 1776	Tube 12X2mm, 5m long (not shown)	1
20	407 0144	O - RING	2

RINSING




figure 6

Parts Listing figure 6

Item	Part No.	Description	Quantity
1	927 1406	Rinse Pump Assy including Tubes Developer/Fixer. (Drain > Rinse Pump and Rinse Pump > Basin), Spouts and Connector.....	1
2	927 0256	Rinse Tank Fixer.	1
3	927 0266	Rinse Tank Developer	1
4	927 0226	Drain Developer (including Gasket)	1
5	927 0216	Drain Fixer (including Gasket)	1
6	430 2063	Screw (Phillips Head, self tapping)	6
7	448 0032	Washer	6
8	428 1421	Screw (Counter Sunk) alternative Pos. 13	4
9	927 0283	Fitting	2
10	450 0701	Lock Washer	2
11	485 3601	Hexagon Nut M12X2mm	2
12	927 1861	Bracket (ReedContact Array Developer/Fixer)	2
13	428 3951	Screw (Hexagon Head, M5X20mm)	4

WATER INLET, DRAIN

NOTE !

The WATER INLET PIPES are installed
with holes pointing down !




figure 7

Parts Listing figure 7

Item	Part No.	Description	Quantity
1	927 0806	Water Valve	1
2	428 3121	Hexagon Screw M4X10	2
3	927 0306	Tube reinforced	2
4	927 0316	Tube Developer/Fixer (Overflow > Drain)	2
5	927 0236	Water Inlet Pipe Fixer	1
6	927 0246	Water Inlet Pipe Developer	1
7	927 0841	Fitting	2
8	917 1181	Hexagon Nut	2
9	927 0826	Adapter 3/4 NH (including Gasket PN 927 0811)	1
10	927 0386	Gasket Set (Drain and Service Keyboard)	1
11	927 0226	Drain Developer (including Gasket)	1
12	927 0216	Drain Fixer (including Gasket)	1
13	917 0491	Spout	2

CLAMPS




figure 8

Parts Listing figure 8

Item	Part No.	Description	Quantity
1	927 1641	Clamp (Size 19/8).....	4
2	927 1631	Clamp (Size 16/8).....	8
3	927 1621	Clamp (Size 15/8).....	4
4	927 1611	Clamp (Size 9,5/7)	4
5	927 1601	Clamp (Size 8/7).....	2


ELECTRONIC BOX (TRANSFORMER)

figure 9

Parts Listing figure 9

Item	Part No.	Description	Quantity
1	927 1046	Transformer (including Connector)	1
2	427 9340	Screw (Round Pan Head M5X15)	4
3	450 0431	Lock Washer (A5.3)	4
4	448 0171	Washer (A5.3)	4
5	485 3031	Hexagon Nut (M5)	4
6	927 1146	Wire (Ground)	1
7	451 0001	U - Terminal	2


ELECTRONIC BOX (POWER and REMOTE CONNECTOR)

figure 10

Parts Listing figure 10

Item	Part No.	Description	Quantity
1	927 1006	Electronic Box (including Ratingplate PN 927 1181 and Gasket PN 917 1786)	1
2	927 1066	Power Switch Socket (including Bracket PN 927 1341 and Counter Sunk Screw PN 428 3981)	1
3	927 1311	Cover Plate.	1
4	927 1056	Quantity Selector Switch Assy (including Torx Screw PN 428 3991 PN 450 0101 and Remote Control Socket)	1
5	927 1071	Voltage Selector for Fuse 5X20mm and Fuse Holder	1
	927 1081	Voltage Selector for Fuse 6.3X32mm and Fuse Holder	1
6	453 4001	Fuse (slow blow 1.25 A)	1
	453 4011	Fuse (slow blow 2.5 A)	1
7	927 1891	Bracket	1
8	924 6271	Cable (including Terminals)	1
9	927 1901	Cap.	1

ELECTRONIC BOX (SERVICE KEYBOARD, HARNESS, CABLES)




figure 11

Parts Listing figure 11

Item	Part No.	Description	Quantity
1	927 1006	Electronic Box (including Ratingplate PN 927 1181 and Gasket PN 917 1786)	1
2	917 2216	Service Keyboard (without Cable, Connector and Gasket) . . .	1
3	428 3011	Plastic Screw	4
4	927 1126	Flat Cable 14 pole	1
5	927 1881	Bracket	1
6	927 1871	Socket (2X) for Reed Contact Array (including Cable)	1
7	430 1951	Torx Screw	4
8	927 1096	16 pin Socket (including Harness and Connector Socket)	1
9	430 2063	Screw (Phillips Head, self tapping)	4
10	927 1271	Spacer	3
11	428 3291	Screw	3


ELECTRONIC BOX (PRINTED CIRCUIT BOARDS)

figure 12

Parts Listing figure 12

Item	Part No.	Description	Quantity
1	927 1006	Electronic Box (including Ratingplate PN 927 1181 and Gasket PN 917 1786)	1
2	927 1321	Insulation Plate.....	1
3	927 1800	PCB A2 (Control).....	1
4	927 1820	PCB A1 (Microcontroller).....	1
5	927 1816	EEPROM Program (Vesion 1.X).....	1
6	927 0996	Cover.....	1

OPERATOR DISPLAY




figure 13

Parts Listing figure 13

Item	Part No.	Description	Quantity
1	428 3971	Screw (Counter Sunk)	4
2	927 0206	Display (Operator Assy including Gasket)	1
3	927 1010	PCB Display	1
4	927 1331	Insulation Plate	1
5	448 5810	Lock Washer	3
6	419 0011	Hexagon Spacer	3
7	485 3111	Self Locking Nut	3
8	927 1116	Flat Cable 10 pole	1

BOARDS, CONNECTORS, CABLES




figure 14

Parts Listing figure 14

Item	Part No.	Description	Quantity
1	927 1046	Transformer (see figure 9)	1
2	927 1096	16 pin Socket with Harness (see figure 11)	1
3	927 1106	16 pole Plug with Harness	1
4	927 1856	Reed Contact Array (including Float, Hexagon Nuts, Cable and Connector)	2
5	927 2216	Service Keyboard (see figure 11)	1
6	927 1126	Flat Cable (14 pole see figure 11)	1
7	927 0206	Operator Display (see figure 13)	1
8	927 1116	Flat Cable (10 pole see figure 13)	1
9	921 1136	Fan Assy (see figure 1)	1
10	927 1056	Switch Socket Assy (see figure 10)	1
11	927 1871	Socket (2X) for Reed Contact Array (including Cable)	1

Part No.	Description	Figure No.
??? ???? ?	GISD Chemistry (not available yet)	
407 0144	O-Ring	5
419 0011	Hexagon Spacer	13
427 9340	Screw (Round Pan Head M5X15)	9
428 1421	Screw (Counter Sunk) alternative Pos. 13	6
428 1461	Screw	1
428 3011	Plastic Screw	11
428 3121	Screw M4X10 Hexagon	7
428 3291	Screw	11
428 3951	Screw	6
428 3961	Screw M6X40 Hexagon	2
428 3971	Screw (Counter Sunk)	13
430 1951	Screw (Torx)	11
430 2063	Screw (Phillips Head, self tapping)	1,6,11
430 2081	Screw	1
430 2101	Screw	1
448 0032	Washer	1,6
448 0171	Washer (A5.3)	9
448 0181	Washer	1
448 0271	Washer	1
448 5810	Lock Washer	13
450 0431	Lock Washer (A5.3)	9
450 0681	Lock Washer	2
450 0701	Lock Washer	6
451 0001	U - Terminal	9
453 4001	Fuse (slow blow 1.25 A)	10
453 4011	Fuse (slow blow 2.5 A)	10
485 2771	Hexagon Nut	1
485 3031	Hex Nut (M5)	9
485 3111	Hexagon Nut (Self Locking)	13
485 3481	Hexagon Nut M6	2
485 3601	Hexagon Nut M12X2mm	6
715 0609	Developer + Fixer Basin Assy RA Chemistry (Gallons)	3
715 7068	Developer + Fixer Basin Assy RP Chemistry (litres)	3
715 7076	Developer + Fixer Basin Assy RP Chemistry (Gallons)	3
917 0041	Rivet Plastic	3
917 0171	O-Ring	2,5

Part No.	Description	Figure No.
917 0266	Tube Developer/Fixer (Mixing Valve >> Floating Lid)	5
917 0276	Tube Developer (Mixing Tank >> Solenoid)	5
917 0306	Tube Fixer (Mixing Tank >> Solenoid)	5
917 0421	Locking Caster	1
917 0431	Caster	1
917 0491	Spout	7
917 0624	Transfer Valve Developer	5
917 0634	Transfer Valve Fixer	5
917 0654	Developer Mixing Valve Assy	5
917 0664	Fixer Mixing Valve Assy	5
917 0681	Valve Bracket	1
917 0704	Bracket	1
917 0824	Elbow Fitting	5
917 0861	Elbow	5
917 1101	Support	1
917 1136	Floating Lid	5
917 1181	Hexagon Nut	7
917 1381	Cover (Transfer Valves)	1
917 1451	Magnet	1
917 2216	Service Keyboard (without Cable Connector and Gasket)	11,14
917 2471	Screw (Phillips)	1
917 6141	Cover (Circulation Pumps)	1
917 6206	Tube Developer/Fixer (Holding Tank >> Output)	5
921 1136	Fan Assy (see figure 1)	14
924 6271	Cable (including Terminals)	10
927 0041	T - Fitting	5
927 0052	Exhaust Tube Assy without O-Ring	2
927 0061	Bracket	2
927 0072	Bottom Plate	2
927 0082	Hexagon Nut (Plastic)	2
927 0096	Gasket	2
927 0106	Tank Developer	1
927 0107	Tank Fixer	1
927 0116	Shroud	1
927 0126	Lower Body Fixer	1
927 0136	Lower Body Developer	1
927 0156	Front Panel Assy with Magnet and Screw	1
927 0166	Rear Panel	1

Part No.	Description	Figure No.
927 0176	Cover	1
927 0191	Shaft	1
927 0206	Display (Operator Assy including Gasket)	13,14
927 0216	Drain Fixer (including Gasket)	6,7
927 0226	Drain Developer (including Gasket)	6,7
927 0236	Water Pipe Fixer	7
927 0246	Water Pipe Developer	7
927 0256	Rinse Tank Fixer	6
927 0266	Rinse Tank Developer	6
927 0281	Screw Cap Rinse	3
927 0282	Coupling Rinse	3
927 0283	Fitting	6
927 0306	Tube reinforced	7
927 0316	Tube Developer/Fixer (Overflow >> Drain)	7
927 0331	Tube	4
927 0386	Gasket Set (Drain and Service Keyboard)	7
927 0406	Mixing Basin Developer	3
927 0416	Mixing Basin Fixer	3
927 0456	Fixer Basin without Accessories	3
927 0466	Developer Basin without Accessories	3
927 0496	Circulation Pump	5
927 0506	Piercer Assy complete with Hexagon Plastic Screw, and Tube Nozzle	3
927 0551	Screw Hexagon Plastic	4
927 0561	Tube Nozzle	4
927 0581	T-Fitting	4
927 0706	Template RP Fixer (litres)	3
927 0716	Template RP Developer (litres)	3
927 0726	Template RP Fixer (Gallons)	3
927 0736	Template RP Developer (Gallons)	3
927 0746	Template RA Fixer (Gallons)	3
927 0756	Template RA Developer (Gallons)	3
927 0806	Water Valve	7
927 0826	Adapter 3/4 NH (including Gasket PN 927 0811)	7
927 0841	Fitting	7
927 0996	Cover	12
927 1006	Electronic Box (including Ratingplate PN 927 1181 and Gasket PN 917 1786)	10,11,12
927 1010	PCB Display	13
927 1046	Transformer (including Connector)	9,14

Part No.	Description	Figure No.
927 1056	Quantity Selector Switch Assy	10,14
927 1066	Power Switch Socket (including Bracket and Conter Sunk Screw)	10
927 1071	Voltage Selector for Fuse 5X20mm and Fuse Holder	10
927 1081	Voltage Selector for Fuse 6.3X32mm and Fuse Holder	10
927 1096	16 pin Socket (including Harness and Connector Socket)	11,14
927 1106	16 pole Plug with Harness	14
927 1116	Flat Cable 10 pole	13,14
927 1126	Flat Cable 14 pole	11,14
927 1136	Fan Assy including Connector	1,14
927 1146	Wire (Ground)	9
927 1271	Spacer	11
927 1311	Cover Plate	10
927 1321	Insulation Plate	12
927 1331	Insulation Plate	13
927 1406	Rinse Pump Assy Developer/Fixer including Tubes	6
927 1601	Clamp (Size 8/7)	8
927 1611	Clamp (Size 9,5/7)	8
927 1621	Clamp (Size 15/8)	8
927 1631	Clamp (Size 16/8)	8
927 1641	Clamp (Size 19/8)	8
927 1701	Tube	4
927 1711	Tube	4
927 1721	Tube	4
927 1736	Tube Fixer (Solenoid >> Circulation Pump)	5
927 1746	Tube Developer (Solenoid >> Circulation Pump)	5
927 1756	Tube Fixer (Circulation Pump >> Mixing Valve)	5
927 1766	Tube Developer (Circulation Pump >> Mixing Valve)	5
927 1776	Tube 12X2mm 5m long (not shown)	5
927 1800	PCB A2 (Control)	12
927 1816	EEPROM Program (Vesion 1.0)	12
927 1820	PCB A1 (Microcontroller)	12
927 1856	Reed Contact Array with Float, Hexagon Nuts, Cable and Connector	1,14
927 1861	Bracket (ReedContact Array)	6
927 1871	Socket (2X) for Reed Contact Array (including Cable)	11,14
927 1881	Bracket	11
927 1891	Bracket	10
927 1901	Cap	10

Description	Part No.	Figure No.
16 pin Socket (including Harness and Connector Socket)	927 1096	11,14
16 pole Plug with Harness	927 1106	14
Adapter 3/4 NH (including Gasket PN 927 0811)	927 0826	7
Bottom Plate	927 0072	2
Bracket	917 0704	1
Bracket	927 0061	2
Bracket	927 1881	11
Bracket	927 1891	10
Bracket (ReedContact Array)	927 1861	6
Cable (including Terminals)	924 6271	10
Cap	927 1901	10
Caster	917 0431	1
Circulation Pump	927 0496	5
Clamp (Size 15/8)	927 1621	8
Clamp (Size 16/8)	927 1631	8
Clamp (Size 19/8)	927 1641	8
Clamp (Size 8/7)	927 1601	8
Clamp (Size 9,5/7)	927 1611	8
Coupling Rinse	927 0282	3
Cover	927 0176	1
Cover	927 0996	12
Cover (Circulation Pumps)	917 6141	1
Cover (Transfer Valves)	917 1381	1
Cover Plate	927 1311	10
Developer + Fixer Basin Assy RA Chemistry (Gallons)	715 0609	3
Developer + Fixer Basin Assy RP Chemistry (Gallons)	715 7076	3
Developer + Fixer Basin Assy RP Chemistry (litres)	715 7068	3
Developer Basin without Accessories	927 0466	3
Developer Mixing Valve Assy	917 0654	5
Display (Operator Assy including Gasket)	927 0206	13,14
Drain Developer (including Gasket)	927 0226	6,7
Drain Fixer (including Gasket)	927 0216	6,7
EEPROM Program (Vesion 1.0)	927 1816	12
Elbow	917 0861	5
Elbow Fitting	917 0824	5
Electronic Box (including Ratingplate PN 927 1181 and Gasket PN 917 1786)	927 1006	10,11,12
Exhaust Tube Assy without O-Ring	927 0052	2

Description	Part No.	Figure No.
Fan Assy (see figure 1)	921 1136	14
Fan Assy including Connector	927 1136	1,14
Fitting	927 0283	6
Fitting	927 0841	7
Fixer Basin without Accessories	927 0456	3
Fixer Mixing Valve Assy	917 0664	5
Flat Cable 10 pole	927 1116	13,14
Flat Cable 14 pole	927 1126	11,14
Floating Lid	917 1136	5
Front Panel Assy with Magnet and Screw	927 0156	1
Fuse (slow blow 1.25 A)	453 4001	10
Fuse (slow blow 2.5 A)	453 4011	10
Gasket	927 0096	2
Gasket Set (Drain and Service Keyboard)	927 0386	7
GISD Chemistry (not available yet)	??? ????	
Hex Nut (M5)	485 3031	9
Hexagon Nut	485 2771	1
Hexagon Nut	917 1181	7
Hexagon Nut (Plastic)	927 0082	2
Hexagon Nut (Self Locking)	485 3111	13
Hexagon Nut M12X2mm	485 3601	6
Hexagon Nut M6	485 3481	2
Hexagon Spacer	419 0011	13
Insulation Plate	927 1321	12
Insulation Plate	927 1331	13
Lock Washer	448 5810	13
Lock Washer	450 0681	2
Lock Washer	450 0701	6
Lock Washer (A5.3)	450 0431	9
Locking Caster	917 0421	1
Lower Body Developer	927 0136	1
Lower Body Fixer	927 0126	1
Magnet	917 1451	1
Mixing Basin Developer	927 0406	3
Mixing Basin Fixer	927 0416	3
O-Ring	407 0144	5
O-Ring	917 0171	2,5
PCB A1 (Microcontroller)	927 1820	12

Description	Part No.	Figure No.
PCB A2 (Control)	927 1800	12
PCB Display	927 1010	13
Piercer Assy complete with Hexagon Plastic Screw, and Tube Nozzle	927 0506	3
Plastic Screw	428 3011	11
Power Switch Socket (including Bracket and Conter Sunk Screw)	927 1066	10
Quantity Selector Switch Assy	927 1056	10,14
Rear Panel	927 0166	1
Reed Contact Array with Float, Hexagon Nuts, Cable and Connector	927 1856	1,14
Rinse Pump Assy Developer/Fixer including Tubes	927 1406	6
Rinse Tank Developer	927 0266	6
Rinse Tank Fixer	927 0256	6
Rivet Plastic	917 0041	3
Screw	428 1461	1
Screw	428 3291	11
Screw	428 3951	6
Screw	430 2081	1
Screw	430 2101	1
Screw (Counter Sunk)	428 3971	13
Screw (Counter Sunk) alternative Pos. 13	428 1421	6
Screw (Phillips Head, self tapping)	430 2063	1,6,11
Screw (Phillips)	917 2471	1
Screw (Round Pan Head M5X15)	427 9340	9
Screw (Torx)	430 1951	11
Screw Cap Rinse	927 0281	3
Screw Hexagon Plastic	927 0551	4
Screw M4X10 Hexagon	428 3121	7
Screw M6X40 Hexagon	428 3961	2
Service Keyboard (without Cable Connector and Gasket)	917 2216	11,14
Shaft	927 0191	1
Shroud	927 0116	1
Socket (2X) for Reed Contact Array (including Cable)	927 1871	11,14
Spacer	927 1271	11
Spout	917 0491	7
Support	917 1101	1
T - Fitting	927 0041	5
T-Fitting	927 0581	4
Tank Developer	927 0106	1
Tank Fixer	927 0107	1

Description	Part No.	Figure No.
Template RA Developer (Gallons)	927 0756	3
Template RA Fixer (Gallons)	927 0746	3
Template RP Developer (Gallons)	927 0736	3
Template RP Developer (litres)	927 0716	3
Template RP Fixer (Gallons)	927 0726	3
Template RP Fixer (litres)	927 0706	3
Transfer Valve Fixer	917 0634	5
Transfer Valve Developer	917 0624	5
Transformer (including Connector)	927 1046	9,14
Tube	927 0331	4
Tube	927 1701	4
Tube	927 1711	4
Tube	927 1721	4
Tube 12X2mm 5m long (not shown)	927 1776	5
Tube Developer (Circulation Pump >> Mixing Valve)	927 1766	5
Tube Developer (Mixing Tank >> Solenoid)	917 0276	5
Tube Developer (Solenoid >> Circulation Pump)	927 1746	5
Tube Developer/Fixer (Holding Tank >> Output)	917 6206	5
Tube Developer/Fixer (Mixing Valve >> Floating Lid)	917 0266	5
Tube Developer/Fixer (Overflow >> Drain)	927 0316	7
Tube Fixer (Circulation Pump >> Mixing Valve)	927 1756	5
Tube Fixer (Mixing Tank >> Solenoid)	917 0306	5
Tube Fixer (Solenoid >> Circulation Pump)	927 1736	5
Tube Nozzle	927 0561	4
Tube reinforced	927 0306	7
U - Terminal	451 0001	9
Valve Bracket	917 0681	1
Voltage Selector for Fuse 5X20mm and Fuse Holder	927 1071	10
Voltage Selector for Fuse 6.3X32mm and Fuse Holder	927 1081	10
Washer	448 0032	1,6
Washer	448 0181	1
Washer	448 0271	1
Washer (A5.3)	448 0171	9
Water Pipe Developer	927 0246	7
Water Pipe Fixer	927 0236	7
Water Valve	927 0806	7
Wire (Ground)	927 1146	9



CUSTOMER EQUIPMENT SERVICES KODAK AG STUTTGART GERMANY